

血细胞分析仪用校准物参考值表

型号：C3

批号：2507C

生产日期：2025-8-11

有效期至：2025-9-24

仪器	参数	单位	靶值	允许偏差
RT-7200	WBC	$\times 10^9/L$	10.0	± 0.3
	RBC	$\times 10^{12}/L$	4.55	± 0.1
	HGB	g/L	147	± 2
	MCV	fL	92.1	± 2
	PLT	$\times 10^9/L$	223	± 15
RT-7300	WBC	$\times 10^9/L$	10.4	± 0.3
	RBC	$\times 10^{12}/L$	4.61	± 0.1
	HGB	g/L	148	± 2
	MCV	fL	88.8	± 2
	PLT	$\times 10^9/L$	249	± 15
RT-7600 RT-7600S	WBC	$\times 10^9/L$	9.9	± 0.3
	RBC	$\times 10^{12}/L$	4.46	± 0.1
	HGB	g/L	148	± 2
	MCV	fL	88.9	± 2
	PLT	$\times 10^9/L$	271	± 15
RT-7900	WBC	$\times 10^9/L$	11.0	± 0.3
	RBC	$\times 10^{12}/L$	4.62	± 0.1
	HGB	g/L	144	± 2
	MCV	fL	84.2	± 2
	PLT	$\times 10^9/L$	257	± 15
RT-7600Vet	WBC	$\times 10^9/L$	10.1	± 0.3
	RBC	$\times 10^{12}/L$	4.25	± 0.1
	HGB	g/L	121	± 2
	MCV	fL	87.5	± 2
	PLT	$\times 10^9/L$	181	± 15