

血细胞分析仪用校准物参考值表

型号：C3

批号：2509C

生产日期：2025-09-30

有效期至：2025-11-13

仪器	参数	单位	靶值	允许偏差
RT-7200	WBC	$\times 10^9/L$	10.16	± 0.3
	RBC	$\times 10^{12}/L$	4.48	± 0.1
	HGB	g/L	135	± 2
	MCV	fL	86.5	± 2
	PLT	$\times 10^9/L$	220	± 15
RT-7300	WBC	$\times 10^9/L$	10.38	± 0.3
	RBC	$\times 10^{12}/L$	4.40	± 0.1
	HGB	g/L	135	± 2
	MCV	fL	83.0	± 2
	PLT	$\times 10^9/L$	235	± 15
RT-7600 RT-7600S	WBC	$\times 10^9/L$	9.93	± 0.3
	RBC	$\times 10^{12}/L$	4.33	± 0.1
	HGB	g/L	136	± 2
	MCV	fL	83.3	± 2
	PLT	$\times 10^9/L$	253	± 15
RT-7900	WBC	$\times 10^9/L$	10.49	± 0.3
	RBC	$\times 10^{12}/L$	4.47	± 0.1
	HGB	g/L	137	± 2
	MCV	fL	79.2	± 2
	PLT	$\times 10^9/L$	251	± 15
RT-7600Vet	WBC	$\times 10^9/L$	10.23	± 0.3
	RBC	$\times 10^{12}/L$	4.19	± 0.1
	HGB	g/L	111	± 2
	MCV	fL	80.8	± 2
	PLT	$\times 10^9/L$	178	± 15